

QDOT MICRO™

Catheter

The Value of QDOT MICRO™ Catheter
Technology



ADVANCEMENT IN RADIOFREQUENCY CATHETER ABLATION OF ATRIAL FIBRILLATION

Radiofrequency (RF) catheters with irrigation and contact force technology are the most advanced and most commonly used technologies for the ablation of atrial fibrillation (AF), providing superior performance as compared to other available modalities.¹⁻⁵ However, technological advancement that enables shortening the duration of ablation via use of high-power (>40W) and short duration (4s) ablations has the potential to further enhance procedural efficiency and improve patient outcomes.⁶

Effective RF ablation of AF is dependent on catheter stability, contact force, power output, ablation time, and temperature.⁴ Measurements of the temperature at the catheter-tissue interface may fluctuate and require operators to carefully control RF power output to ensure safe tissue temperatures and achieve optimal outcomes.⁵



Standard RF ablation procedures are conducted at 20W-40W,^{7,8} which may lead to lengthy **procedure times**.¹⁻³



Lengthy application of low RF power may **increase the risk of damage to surrounding tissues**, such as the esophagus, which may increase the risk of atrioesophageal fistula.⁷⁻¹¹



Fluctuations in temperature measurement and conservative use of RF power may **increase the risk of imprecise lesion formation**.^{7,12-14}


QDOT MICRO™ TECHNOLOGY: INNOVATION IMPROVING OUTCOMES

QDOT MICRO™ Technology enables automation of key ablation parameters to standardize and optimize RF energy delivery, which is enabled by thermocouples providing real-time feedback on catheter surface temperature. This technology provides temperature control, enables a high power, short duration workflow and allows for high-resolution substrate mapping capabilities.

OPTIMIZED TEMPERATURE CONTROL

QDOT MICRO™ Technology has been designed to **provide temperature control through automatic adjustment of power and fluid output** based on real-time temperature measurement. Use of higher average RF power with QDOT MICRO™ Technology is safely enabled by this precise temperature control.^{1-7;15-16}

HIGH POWER, SHORT DURATION ABLATION



The temperature control of QDOT MICRO™ Technology allows for the **safe use of higher RF power (up to 90 Watts)** in short bursts (up to 4 seconds) in QMODE+™ Ablation Mode to improve ablation efficiency without compromising safety.^{7;11;16-18}

Compared to conventional contact force sensing catheters, high-power short duration (vHP-SD, 90 W/4 s) QDOT MICRO™ Catheter (QMODE+™ Ablation Mode based PVI) provides comparably safe and effective PVI.^{6;15-18}

VALUE OF QDOT MICRO™ TECHNOLOGY

Innovative QDOT MICRO™ Technology enables improved procedure efficiency, allows for very low complication rates, facilitates the reduction of fluid delivery and fluoroscopy time, while providing excellent efficacy.

IMPROVED EFFICIENCY




60 MIN
PROCEDURE DURATION
WITH QMODE+™
ABLATION MODE^{6;16-20}

Ablation with the QDOT MICRO™ Catheter has demonstrated a **median procedure duration of ~60 minutes** as compared to a median ~105 minutes with conventional CF-sensing catheters.^{6;16-20}



UP TO 49%
SHORTER PROCEDURE
TIME^{2-3;11;16;21-22}

Ablation using QDOT MICRO™ Catheter can **reduce procedure time by up to 49%** as compared to standard contact force catheters.^{1-3;11;16;23-24*} QDOT MICRO™ Technology is also compatible with VISITAG SURPOINT™ Module (i.e. Ablation Index), which may allow for a greater potential reduction in procedure time.^{6;22}



Comparison of QDOT MICRO™ Technology to other studies of conventional contact force (CF) catheters for total procedure time.

* Data for conventional CF catheters is the average (204.7 min) of THERMOCOOL SMARTTOUCH® SMART-AF (Natale *et al*, 2014; n=160, 222.0 min), THERMOCOOL SMARTTOUCH® SMART-SF (Chinitz *et al*, 2018; n=159, 181.1 min), and TactiCath TOCCATA (Reddy *et al*, 2012; n=34, 211.0 min).

** QDOT QEFFICIENCY (Osorio *et al* 2022) reported a median procedure time, whereas all other studies reported mean procedure times.

† QDOT QEFFICIENCY (Osorio *et al*, 2022), and TactiCath TOCCASTAR (Reddy *et al*, 2015), studies report a median procedure time, whereas the other studies reported in this graph report mean procedure times.^{6;11;16;22;29}


REDUCED ABLATION TIME



**UP TO 62%
SHORTER ABLATION
TIME WITH QMODE+™
ABLATION MODE**
11;16;21-22;24-25

Total ablation time with QDOT MICRO™ Technology in QMODE+™ (very high power, short duration RF bursts) was **up to 62% shorter** than with conventional contact force catheter ablation technologies^{11;16;21-22;24-25}

† Based on a weighted average of the procedure time with QDOT MICRO™ Technology in QMODE+™ Ablation Mode (129.8 min; N = 42) and in QMODE+™ Ablation Mode (105.2 min; N = 52), and the weighted average of the procedure time with irrigated, contact force RF catheters (202.8 min; 5 studies, N = 622).^{1-3;11;26-27}
The third-party trademarks used herein are the trademarks of the respective owners.



Comparison of QDOT MICRO™ Technology to other studies of conventional contact force (CF) catheters for total ablation time.

* Data for conventional CF catheters is the average (122.6 min) of THERMOCOOL SMARTTOUCH® SMART-AF (Natale *et al.*, 2014; n=160, 121.5 min), THERMOCOOL SMARTTOUCH® SMART-SF (Chinitz *et al.*, 2018; n=159, 104.3 min), TactiCath TOCCATA (Reddy *et al.*, 2012; n=34, 218.0 min), and TactiCath TOCCASTAR (Reddy *et al.*, 2015; n=152, 46.5 min).


** QDOT QEFFICIENCY (Osorio *et al.*, 2022) and TactiCath TOCCASTAR (Reddy *et al.*, 2015) studies reported a median total ablation time, whereas all other studies reported mean total ablation times.^{5;11;16;22;29}

REDUCED RF APPLICATION TIME



**UP TO 84%
REDUCTION IN RF
APPLICATION TIME
WITH QMODE+™
ABLATION MODE**
3;11;16;21-22

In the QDOT QEFFICIENCY trial, QDOT Technology had **up to 84% shorter RF application time** than conventional contact force catheters.^{3;11;16;21;28}



Comparison of QDOT MICRO™ Technology to other studies of conventional contact force (CF) catheters for RF application time.

* Data for conventional CF catheters is the average (49.2 min) of THERMOCOOL SMARTTOUCH® SMART-AF (Reddy *et al.*, 2019, n=160, 60.6 min), THERMOCOOL SMARTTOUCH® SMART-SF (Chinitz *et al.*, 2018; n=159, 49.5 min), and TactiCath TOCCATA (Reddy *et al.*, 2012, n=34, 37.6 min).


** QDOT QEFFICIENCY (Osorio *et al.*, 2022)) reported a median RF application time, whereas all other studies reported mean RF application times.^{5;11;16;22;29}

REDUCED FLUOROSCOPY EXPOSURE



**UP TO 79%
SHORTER MEAN
FLUOROSCOPY TIME**
1-3;11;16;19;21-22;27;29;33

RF ablation with QDOT MICRO™ technologies had **up to 79% shorter mean fluoroscopy time** compared to conventional contact force catheters.^{1-3;11;16;19;21-22;26-27}



Comparison of QDOT MICRO™ Technology to other studies of conventional contact force (CF) catheters for fluoroscopy time.

* Data for conventional CF catheters is the average (30.9 min) of THERMOCOOL SMARTTOUCH® SMART-AF (Natale *et al*, 2014; n=160, 41.5 min), THERMOCOOL SMARTTOUCH® SMART-SF (Chinitz *et al*, 2018; n=159, 18.6 min), TactiCath TOCCATA (Reddy *et al* 2012; n=34, 36.3 min), and TactiCath TOCCASTAR (Reddy *et al*, 2015, n=150, 27.0 min).

** QDOT QEFFICIENCY (Osorio *et al*, 2022) and TactiCath TOCCASTAR (Reddy *et al*, 2015) studies reported a median fluoroscopy time, whereas all other studies reported mean fluoroscopy times.^{6;11;16;22;29}



**PROCEDURES
PERFORMED UNDER
MILD CONSCIOUS
SEDATION¹⁹**

In a retrospective study 90% of patients treated with very high-power short duration (vHPSD) ablations were **successfully performed under mild conscious sedation.**¹⁹

LOW COMPLICATION RATE

Very high-power, short duration ablation **may help to reduce the incidence of complications** associated with poor catheter stability and tip-to tissue contact.¹¹



**EXCELLENT SAFETY
PROFILE WITH QDOT
MICRO™ TECHNOLOGY¹⁶⁻¹⁹**

In prospective, non-randomized, multicenter studies, **catheter ablation with QDOT MICRO™ Technology** was associated with no incidence of steam pops or charring. Use of QMODE+™ Ablation Mode was associated with **no incidence of stroke, atriopharyngeal fistula, and PV stenosis.**^{1;2;6;11;23;26}

REDUCED FLUID DELIVERY


Excessive fluid delivery during RF ablation can increase the risk of complications, including fluid overload, heart failure, acute exacerbation of heart failure, acute respiratory distress, hypoxia, and pulmonary edema.³⁵⁻³⁸



**UP TO 72%
LESS FLUID DELIVERED
WITH QMODE MICRO™
TECHNOLOGY**^{11;16-17;21-22;29}

QDOT MICRO™ Technology delivered 72% less fluid with QMODE+™ Ablation Mode and **59% less fluid** with QMODE™ Ablation Mode than conventional contact-force catheter ablation technology.^{11;16-17;21-22;26*†}

Delivery of less fluid during ablation may reduce the risk of complications, particularly in patients with a high comorbidity burden.^{21†}



Comparison of QDOT MICRO™ Technology to other studies of conventional contact force (CF) catheters for fluid volume.

* Data for conventional CF catheters is the average (1389.0 mL) of THERMOCOOL SMARTTOUCH® SMART-AF (Natale *et al*, 2014; n=158, 1879.6 mL) and THERMOCOOL SMARTTOUCH® SMART-SF (Chinitz *et al*, 2018; n=156, 898.4 mL).

** QDOT QEFFICIENCY (Osorio *et al*, 2022) reported a median fluid volume, whereas all other studies reported mean fluid volumes.^{6;11;16;22;29}

EFFICACY



**HIGHER
PROBABILITY OF FIRST
PASS ISOLATION**³⁹

A prospective study (156 patients) using very high-power short duration (vHPSD) ablation was associated with a higher probability of first-pass isolation compared to low-power long-duration (LPLD) ablation (OR = 2.90, $p = 0.014$).³⁹

FREEDOM FROM ARRHYTHMIA AT 12 MONTHS



**86%
OF PATIENTS ARE
ARRHYTHMIA-FREE AT
12 MONTHS**¹⁶

In a retrospective analysis, **86.6% of patients were free from recurrence of atrial arrhythmia** after ablation procedures using QMODE+™ ablation mode.^{11;16;18-20;29}

REFERENCES

1. Chinitz LA, Melby DP, Marchlinski FE, Delaughter C, Fishel RS, *et al.* (2017) Safety and efficiency of porous-tip contact-force catheter for drug-refractory symptomatic paroxysmal atrial fibrillation ablation: results from the SMART SF trial. *J Ep Europace* 20 (FI_3): f392-f400.
2. Natale A, Reddy VY, Monir G, Wilber DJ, Lindsay BD, *et al.* (2014) Paroxysmal AF catheter ablation with a contact force sensing catheter: results of the prospective, multicenter SMART-AF trial. *J Am Coll Cardiol* 64 (7): 647-656.
3. Reddy VY, Shah D, Kautzner J, Schmidt B, Saoudi N, *et al.* (2012) The relationship between contact force and clinical outcome during radiofrequency catheter ablation of atrial fibrillation in the TOCCATA study. *Heart Rhythm* 9 (11): 1789-1795.
4. Calkins H, Hindricks G, Cappato R, Kim Y-H, Saad EB, *et al.* (2017) 2017 HRS/EHRA/ECAS/APHRS/SOLAECE expert consensus statement on catheter and surgical ablation of atrial fibrillation. *J Ep Europace* 20 (1): e1-e160.
5. Stevenson WG, Cooper J, Sapp J (2004) Optimizing RF output for cooled RF ablation. *J Journal of cardiovascular electrophysiology* 15 S24-S27.
6. Tilz R.R., Sano M., Vogler J., Fink T., *et al.* (2021). Very high-power short-duration temperature-controlled ablation versus conventional power-controlled ablation for pulmonary vein isolation: The fast and furious - AF study. *IJC Heart & Vasculature*. Vol. 35: 100847
7. Leshem E, Zilberman I, Tschabrunn CM, Barkagan M, Contreras-Valdes FM, *et al.* (2018) High-power and short-duration ablation for pulmonary vein isolation: biophysical characterization. *4 (4): 467-479.*
8. Barkagan M, Contreras-Valdes FM, Leshem E, Buxton AE, Nakagawa H, *et al.* (2018) High-power and short-duration ablation for pulmonary vein isolation: Safety, efficacy, and long-term durability. *Journal of Cardiovascular Electrophysiology* 29 (9): 1287-1296.
9. Leshem E, Tschabrunn CM, Contreras-Valdes FM, Zilberman I, Anter E (2017) Evaluation of ablation catheter technology: comparison between thigh preparation model and an in vivo beating heart. *Heart Rhythm* 14 (8): 1234-1240.
10. Haines DE (1993) The biophysics of radiofrequency catheter ablation in the heart: the importance of temperature monitoring. *Pacing and Clinical Electrophysiology* 16 (3): 586-591.
11. Reddy VY, Grimaldi M, De Potter T, Vijgen JM, Bulava A, *et al.* (2019) Pulmonary Vein Isolation with Very High Power-Short Duration Temperature-Controlled Lesions: The First-in-Human QDOT-FAST Multicenter Trial. *JACC: Clinical Electrophysiology* 912.
12. Anter E, Contreras-Valdes FM, Shvilkin A, Tschabrunn CM, Josephson ME (2014) Acute pulmonary vein reconnection is a predictor of atrial fibrillation recurrence following pulmonary vein isolation. *Journal of Interventional Cardiac Electrophysiology* 39 (3): 225-232.
13. Anter E, Tschabrunn CM, Contreras-Valdes FM, Buxton AE, Josephson ME (2014) Radiofrequency ablation annotation algorithm reduces the incidence of linear gaps and reconnection after pulmonary vein isolation. *Heart Rhythm* 11 (5): 783-790.
14. Andrade JG, Monir G, Pollak SJ, Khairy P, Dubuc M, *et al.* (2014) Pulmonary vein isolation using "contact force" ablation: the effect on dormant conduction and longterm freedom from recurrent atrial fibrillation--a prospective study. *Heart Rhythm* 11 (11): 1919-1924.
15. Stabile G., *et al.* (2021). In vivo biophysical characterization of very high power, short duration, temperature-controlled lesions. *Pacing Clin Electrophysiol*. 1-7
16. Osorio J, Hussein AA, Delaughter MC, *et al.* (2022). Very high-power short-duration, temperature-controlled radiofrequency ablation in paroxysmal atrial fibrillation: the prospective multicenter Q-FICIENCY trial. *JACC Clin Electrophysiol*. In press
17. Almorad A, Wielandts JY, Haddad ME, *et al.* (2021) Performance and Safety of Temperature- and Flow-Controlled Radiofrequency Ablation in Ablation Index-Guided Pulmonary Vein Isolation. *JACC Clin Electrophysiol* 7(3): 408-409.
18. Bortone A, Albenque JP, Ramirez D, *et al.* (2022) 90 vs 50-Watt Radiofrequency Applications for Pulmonary Vein Isolation: Experimental and Clinical Findings. *Circ Arrhythm Electrophysiol*. 15:e010663.
19. Chu G, Calvert P, Sidhu B, *et al.* (2022) Patient experience of very high power short duration radiofrequency ablation for atrial fibrillation under mild conscious sedation. *Journal of Interventional Cardiac Electrophysiology*: 1-9.
20. Salló Z, Perge P, Balogi B, *et al.* (2022) Impact of High-Power and Very High-Power Short-Duration Radiofrequency Ablation on Procedure Characteristics and First-Pass Isolation During Pulmonary Vein Isolation. *Front. Cardiovasc. Med*. 9:935705.
21. Chinitz L, Goldstein LJ, Barnow A, Maccioni S, Daskiran M, *et al.* (2018) Comparison of real-world clinical and economic outcomes between the ThermoCool((R)) SF and ThermoCool((R)) catheters in patients undergoing radiofrequency catheter ablation for atrial fibrillation. *Clinicoecon Outcomes Res* 10 587-599.
22. De Potter T., *et al.* (2021). Temperature-Controlled Catheter Ablation for Paroxysmal Atrial Fibrillation: The QDOT-MICRO Workflow Study *Journal of Atrial Fibrillation*. Vol. 13(6): 1-7.
23. DePotter T, Hunter T, Boo L, *et al.* (2020) The industrialization of ablation: a highly standardized and reproducible workflow for radiofrequency ablation of atrial fibrillation. *J Interv Card Electrophysiol* 59: 21-27.
24. UG-5050-017M (00A) CARTO® 3 System QDOT MICRO™ Module IFU. November 2019
25. Seidl S, Müllleder T, Kaiblinger J, *et al.* (2022) Very High-Power Short-Duration (HPSD) Ablation for Pulmonary Vein Isolation: Short and Long-Term Outcome Data. *J. Cardiovasc. Dev. Dis*. 9; 278.
26. De Potter T, Grimaldi M, Jensen HK, Kautzner J, Neuzil P, *et al.* (2019) Acute safety and performance of a novel catheter with optimized temperature control and microelectrodes in the treatment of patients with paroxysmal atrial fibrillation. *European Hearth Rhythm Association*. Lisbon, March 19, 2019.
27. Taghji P, El Haddad M, Philips T, Wolf M, Knecht S, *et al.* (2018) Evaluation of a strategy aiming to enclose the pulmonary veins with contiguous and optimized radiofrequency lesions in paroxysmal atrial fibrillation: a pilot study. *JACC: Clinical Electrophysiology* 4 (1): 99-108.
28. Yavin HD, Leshem E, Shapira-Daniels A, *et al.* (2020). Impact of High-Power Short-Duration Radiofrequency Ablation on Long-Term Lesion Durability for Atrial Fibrillation Ablation *J Am Coll Cardiol EP*: 6:973-85.
29. De Potter T, Grimaldi M, Jensen HK, Kautzner J, Neuzil P, *et al.* (2019) Results Of The QDOT MICRO Study For Paroxysmal Atrial Fibrillation Treatment. *Heart Rhythm*. San Francisco, May 11, 2019.
30. Ullah W, McLean A, Tayebjee MH, Gupta D, Ginks MR, *et al.* (2016) Randomized trial comparing pulmonary vein isolation using the SmartTouch catheter with or without real-time contact force data. *Heart Rhythm* 13 (9): 1761-1767.
31. Gupta D, De Potter T, Disher T, *et al.* (2020). Comparative effectiveness of catheter ablation devices in the treatment of atrial

fibrillation: a network meta-analysis. *J. Comp. Eff. Res.* 2020 Jan;9(2):115-126. doi: 10.2217/ce-2019-0165.

32. Klein G, Lickfett L, Schreieck J, Deneke T, Wiczorek M, *et al.* (2015) Comparison of 'anatomically designed' and 'point-by-point' catheter ablations for human atrial fibrillation in terms of procedure timing and costs in German hospitals. *Ep Europace* 17 (7): 1030-1037.
33. Reddy VY, Dukkipati SR, Neuzil P, Natale A, Albenque JP, *et al.* (2015) Randomized, Controlled Trial of the Safety and Effectiveness of a Contact Force-Sensing Irrigated Catheter for Ablation of Paroxysmal Atrial Fibrillation: Results of the TactiCath Contact Force Ablation Catheter Study for Atrial Fibrillation (TOCCASTAR) Study. *Circulation* 132 (10): 907-915.
34. Kirchhof P, Benussi S, Kotecha D, Ahlsson A, Atar D, *et al.* (2016) 2016 ESC Guidelines for the management of atrial fibrillation developed in collaboration with EACTS. *J European heart journal* 37 (38): 2893-2962.
35. Hwang ES, Pak HN, Park SW, Park JS, Joung B, *et al.* (2010) Risks and benefits of an open irrigation tip catheter in intensive radiofrequency catheter ablation in patients with non-paroxysmal atrial fibrillation. *Circ J* 74 (4): 644-649.
36. Waldo AL, Wilber DJ, Marchlinski FE, Stevenson WG, Aker B, *et al.* (2012) Safety of the open-irrigated ablation catheter for radiofrequency ablation: safety analysis from six clinical studies. *Pacing Clin Electrophysiol* 35 (9): 1081-1089.
37. Oza SR, Hunter TD, Biviano AB, Dandamudi G, Herweg B, *et al.* (2014) Acute safety of an open-irrigated ablation catheter with 56-hole porous tip for radiofrequency ablation of paroxysmal atrial fibrillation: analysis from 2 observational registry studies. *J Cardiovasc Electrophysiol* 25 (8): 852-858.
38. Huang HD, Waks JW, Contreras-Valdes FM, Haffajee C, Buxton AE, *et al.* (2016) Incidence and risk factors for symptomatic heart failure after catheter ablation of atrial fibrillation and atrial flutter. *Europace* 18 (4): 521-530.
39. Salló Z, Perge P, Balogi B, *et al.* (2022) Impact of High-Power and Very High-Power Short-Duration Radiofrequency Ablation on Procedure Characteristics and First-Pass Isolation During Pulmonary Vein Isolation. *Front. Cardiovasc. Med.* 9:935705.

QDOT MICRO™ Uni-Directional Navigation Catheter Biosense Webster

The Biosense Webster QDOT MICRO™ Uni-Directional Catheter is a steerable multi-electrode luminal catheter with a deflectable tip designed to facilitate electrophysiological mapping of the heart and to transmit radiofrequency (RF) energy to the catheter tip electrode for ablation purposes.

Indications for Use

The Biosense Webster QDOT MICRO™ Uni-Directional Navigation Catheter and related accessory devices are indicated for catheter-based cardiac electrophysiological mapping (stimulating and recording) and, when used in conjunction with a compatible radiofrequency generator, for cardiac ablation. The Biosense Webster QDOT MICRO™ Uni-Directional Navigation Catheter provides a real-time measurement of temperature and contact force between the catheter tip and heart wall, as well as location information when used with the CARTO™ 3 System.

Contraindications

Do not use this catheter:

1. If the patient has had a ventriculotomy or atriotomy within the preceding eight weeks because the recent surgery may increase the risk of perforation.
2. In patients with a myxoma or an intracardiac thrombus as the catheter could precipitate an embolus.
3. In patients with prosthetic valves as the catheter may damage the prosthesis.
4. In the coronary arterial vasculature due to risk of damage to the coronary arterial vasculature.
5. In patients with an active systemic infection because this may increase the risk of cardiac infection.
6. Via the transseptal approach in a patient with an interatrial baffle or patch because the opening could persist and produce an iatrogenic atrial shunt.
7. Via the retrograde trans-aortic approach in patients who have had aortic valve replacement.
8. With a long sheath or short introducer < 8.5 F in order to avoid damage to the catheter shaft.

Warnings and Precautions

Do not use excessive force to advance or withdraw the catheter when resistance is encountered during catheter manipulation through the sheath. Do not manually pre-shape the distal shaft of the catheter by applying external forces intended to bend or affect the intended shape or curve of the catheter. Prior to use, the catheter must be warmed up as specified in the Directions for Use section. If the catheter has not reached a steady state condition, there is potential for a zero-offset drift to occur which could result in an inaccurate contact force reading.

For all product details please consult the IFU (instruction for use) of this product.

For product details such as indications, contraindications, warnings and precautions please consult the IFU. This publication is not intended for distribution outside of the EMEA region.

This summary has been written by Biosense Webster (Europe), a division of Johnson & Johnson Medical NV/SA based on the referenced article, and is provided for information purposes only.

Important information: Prior to use, refer to the instructions for use supplied with the device for indications, contraindications, side effects, warnings and precautions.

Biosense Webster
Johnson & Johnson
Walterovo náměstí 329/1
158 00, Praha 5
Czech republic

For additional medical information request, please contact: <https://www.jnjmedtech.com/en-EMEA/mir>

