

# Medtronic

## PULSED AF Pivotal Trial Overview



PULSED AF evaluated the safety and effectiveness of the PulseSelect™ Pulsed Field Ablation System for the treatment of patients with paroxysmal or persistent atrial fibrillation. PulseSelect™ Pulsed Field Ablation System is an investigational device that has not been approved for commercial use in EU.



0 Esophageal events



0 PV stenosis



0 Phrenic nerve injury



0 Coronary artery spasm



1/300 Stroke

0/300 Pericarditis

1/300 Tamponade

0/300 Vagal nerve injury

0/300 Transient ischemic attack

0/300 Vascular access complications

0/300 Major bleeding

0/300 Cardiovascular Hospitalization

0/300 Myocardial infarction

0/300 Death

### Quality of life

Improved significantly and clinically meaningful



**AFEQT score improved by 29.4 and 29.0** (95% CI, 25.5 to 32.5) points in the paroxysmal and persistent populations respectively from baseline to 12 months.



**EQ-5D-5L score improved by 0.05** points in paroxysmal **and 0.06** points in persistent atrial fibrillation patients.

# Trial design & study population

## Trial design


Prospective controlled paired single arm clinical study

### Global multicenter study

**9 countries:** Australia, Austria, Belgium, Canada, France, Japan, Netherlands, Spain, United States.

**41 sites**

**67 operators**



**150 Paroxysmal**  
**150 Persistent**





The study population included 50% (150) recurrent symptomatic paroxysmal, and 50% (150) persistent atrial fibrillation patients refractory to class I or III anti-arrhythmic drugs.



**96% (287)** of patients completed 12-month follow up.



## Freedom from AF/AT/AFL




# Outcomes



## Rigorous arrhythmia monitoring



## Primary effectiveness



Acute procedure failure, AF/AFL/AT recurrence, cardioversion, repeat ablation, new/re-initiated/increased AADs, any subsequent AF surgery

# Efficiency

## Time between first and last application



**Paroxysmal AF**  
58 ± 28 minutes



**Persistent AF**  
64 ± 28 minutes

**Procedure times under 50 minutes**  
when excluding the 20 minute trial-mandated wait period.



**Total PFA energy delivery under 30 seconds**

**67** operators

**91%** first use of the PulseSelect system

## Medtronic

**Europe**  
Medtronic International Trading Sàrl.  
Route du Molliau 31  
Case postale  
CH-1131 Tolochenaz  
www.medtronic.eu  
Tel: +41 (0)21 802 70 00  
Fax: +41 (0)21 802 79 00

Caution: The products are limited to investigational use only and are not approved for commercial use in any geography. CE Mark pending.

UC202310980-update-global-pulsed-af-brochure-en-we-9010920 © Medtronic 2023.  
Medtronic, Medtronic logo, and Engineering the extraordinary are trademarks of Medtronic.  
All other brands are trademarks of a Medtronic company.  
All rights reserved.