

Setting a  
new standard  
in safety.

**Engineered for efficiency.**

**PulseSelect™**

Pulsed Field Ablation (PFA) System

**Medtronic**

# PulseSelect™

Pulsed Field Ablation System

## Setting a new standard in safety

### Differentiated system safety features

Developed over 15 years of PFA research

Automatic overcurrent detection for safe energy delivery

R-wave gating for synchronisation of energy delivery

Test pulse for proximity detection to phrenic nerve



### Unmatched safety

Engineered with differentiated safety features from 15 years of PFA research, and backed by the results of one of the safest IDE AF ablation trials to date.

### Consistent efficiency

Rapid, effective PVI through consistent and predictable energy delivery and catheter maneuverability.

### Simplified adaptability

Seamless transition to PFA with freedom to adapt to your preferred workflow.

Indicated for use for paroxysmal (PAF) and persistent (PsAF) patients.

## PULSED AF<sup>2</sup>

### Trial design

Paired single arm. Prospective, non-randomized, Global IDE clinical study.

**9 countries**

**41 sites**

**67 operators**

**300 PFA subjects**

# 0.7%

### primary safety event rate

One of the lowest safety event rates of any atrial fibrillation technology to-date.

0 Esophageal events

0 PV stenosis

0 Phrenic nerve injury

0 Coronary spasm noted

Total of 13 adverse events measured, resulted in 1 cerebrovascular accident and 1 tamponade.

# Engineered for efficiency

## PulseSelect™ PFA Catheter

9 electrodes built to **sense, ablate, and pace**


25 mm diameter loop

9 Fr shaft with bidirectional steering

**Fixed electrode spacing** to produce predictable, consistent, and contiguous energy delivery


**20-degree forward tilt** to ensure consistent uniform tissue contact

Full catheter visualisation on existing mapping systems and imaging tools<sup>1</sup>




## Intuitive stepwise approach to PVI†

All pulmonary veins isolated in less than 30 seconds of total energy delivery<sup>1</sup>



**Over-the-wire**  
for enhanced stability and self-centering control



**Ostial applications**  
Stepwise catheter rotation for overlapping lesions



**Antral applications**  
Stepwise catheter rotation for overlapping lesions

## Proven efficacy Freedom from AF/AFL/AT

70%  
PAF

62%  
PsAF

## Consistent efficiency

Left atrial dwell time **50 minutes or under** when excluding the 20-minute trial-mandated wait period.

65  
(± 29) minutes  
PAF

70  
(± 29) minutes  
PsAF

## 1st

PFA IDE trial completed


## Only

PFA IDE trial to demonstrate efficacy for **both PAF and PsAF patients.**



## FlexCath Contour™ Steerable Sheath – 10 Fr

### Steering therapy with confidence



Braided bidirectional shaft design provides catheter support and resistance to kinking

Smooth sheath-dilator transition with transseptal crossing indication

Robust adjustment knob and rounded handle

Two tip lengths enable alignment with RIPV in difficult anatomies

## FlexCath Cross™ Transseptal Solution

### Experience the full potential of a zero exchange workflow

#### Zero sheath exchanges

The integrated needle and dilator design removes the need for needle, guidewire, transseptal, and therapeutic sheath exchanges.

#### Versatility for your procedures

Compatible with the industry-leading introducers and with both mechanical and RF crossing.

#### Cross with control

The needle only deploys when activated, designed for ease of control during transseptal puncture.

<sup>1</sup>Operators targeted pulmonary vein isolation using a wide antral approach during PULSED AF procedures.

<sup>1</sup> Medtronic data on file. November 2023.

<sup>2</sup> Verma A, Haines DE, Boersma LV, et al. Pulsed Field Ablation for the Treatment of Atrial Fibrillation: PULSED AF Pivotal Trial. *Circulation*. May 9, 2023;147(19):1422-1432.

CE0123

# Medtronic

#### Europe

Medtronic International Trading Sàrl.  
Route du Molliau 31  
Case postale  
CH-1131 Tolochenaz  
medtronic.eu  
Tel: +41 0 21 802 70 00  
Fax: +41 0 21 802 79 00

medtronic.eu

#### United Kingdom/Ireland

Medtronic Limited  
Building 9  
Croxley Green Business Park  
Hatters Lane  
Watford  
Herts WD18 8WW  
medtronic.co.uk  
Tel: +44 0 1923 212213  
Fax: +44 0 1923 241004

#### Brief Statement

Product CE mark according to law.

See the device manual for information regarding the instructions for use, indications, contraindications, warnings, precautions, and potential adverse events. For further information, contact your local Medtronic representative and/or consult the Medtronic website at [www.medtronic.eu](http://www.medtronic.eu).

For applicable products, consult instructions for use on [manuals.medtronic.com](http://manuals.medtronic.com). Manuals can be viewed using a current version of any major internet browser. For best results, use Adobe Acrobat® Reader with the browser.

UC202402764a-pulseselect-solution-brochure-en-emea-11429877 ©2023 Medtronic. Adobe and Acrobat Reader are registered trademarks of Adobe Systems incorporated in the United States and/or other countries.™™™Third-party brands are trademarks of their respective owners. All other brands are trademarks of a Medtronic company. All Rights Reserved. Printed in Europe.